CIAT
QK
866
.764
C.4

**POLYACRYLAMIDE GEL ELECTROPHORESIS PROCEDURES
FOR CULTIVAR IDENTIFICATION OF FIELD BEAN,
CASSAVA AND PASTURE LEGUMES**

A. HUSSAIN, W. BUSHUK, H. RAMIREZ, W. M. ROCA

Polyacrylamide gel electrophoresis procedures for cultivar identification of field bean(Phaseolus vulgaris), cassava(Manihot esculenta) and the pasture legumes(Centrosema macrocarpum, C. pubescens, C. acutifolium, Desmodium ovalifolium, Stylosanthes capitata).

A. Hussain, W. Bushuk.

Food Science Department,
University of Manitoba.
Winnipeg, Canada, R3T 2N2.
(Publication No. 115)

and

H. Ramirez, W.M. Roca,

Biotechnology Research Unit,
Centro Internacional de Agricultura Tropical-CIAT,
A.A. 6713, Cali, Colombia. S.A.

INTRODUCTION

Morphological descriptors of plant cultivars often present problems in clearcut identification because morphological differences within species are too ambiguous to discriminate. Ideally the differences between cultivars should be based on direct comparison of genes but this approach is still expensive and time consuming. However these differences can be measured indirectly by comparing the the products of the gene activity i.e. by using proteins as genotype markers.

Electrophoretic techniques outlined in the following pages are based on the differential separation of protein molecules varying in net charge and size through polyacylamide gel matrix under the influence of an applied electric field. Proteins encoded by different genes may have different molecular charge or molecular size and thus may have different mobilities in a given electrical field.

It should be noted that a particular electrophoregram does not reveal all the elements of actual genetic structure. The absence of a particular protein band does not necessarily mean that the genetic information required for its synthesis is also absent. It is possible that such information is blocked or switched off in a particular species or tissue. Nevertheless it is now well established that electrophoretic analyses can contribute useful information where morphological diversities are hard to establish.

Field Bean (Phaseolus vulgaris)

Procedure

Preparation of Sample:

1. Remove seed coat with a scalpel.
2. Grind individual seed into a fine powder using mortar and pestle.
3. Weigh 0.1g of seed meal using sensitive balance to nearest milligram.
4. Add 0.325 ml of extraction solution* and vortex for 2-3 min.
5. Incubate the mixture at 40°C for 2h.
6. Centrifuge at 27,000 Xg for 20 min at 10°C.
7. Remove 100 uL of supernatant into a vial.
8. Add 100 uL of dye solution** and mix.

Preparation of gel:

9. Weigh:

Acrylamide	10 g
Sucrose	5 g
Bis-acrylamide	0.5 g
Ascorbic acid	0.1 g
Ferrous sulphate	0.0015 g

10. Add 95 ml of Aluminum lactate buffer, pH 3.1.***
11. Mix with magnetic stirrer and filter; collect filtrate in a conical flask and degas for 5 min.
12. Stand filtrate in freezer compartment of a fridge for 10 min.
13. Add 175 uL of 3% Hydrogen peroxide solution, swirl few times and pour into Bushuk's apparatus as directed. (CIAT Working Document No. 19, 1986-Spanish; Polyacrylamide gel electrophoresis method for wheat variety identification. I.S.O. standard reference method-English)

- * Extraction solution= 0.1M acetic acid containing 0.1% 2-mercaptoethanol.
- ** Dye solution= 0.5% methyl green dye in 40% sucrose solution.
- *** Aluminum lactate buffer= Dissolve 2.5g Aluminum lactate in 950 ml of distilled water, mix and adjust to pH 3.1 with Lactic acid. Finally adjust volume to 1L with dist. water.

Electrophoresis:

Apparatus:	Bushuk's Vertical Unit
Acrylamide concentration:	10%
Running buffer:	Aluminum lactate, pH 3.1
Running current:	30 mA/gel
Running time:	5 h
Running temperature:	20 °C

Staining:

14. After electrophoresis, stain the gel for at least 8 h or overnight with Coomassie Blue R-250.*
15. Rinse the stained gel with distilled water and then wash excess stain with soapy water, rinse again and then destain with 12% Trichloro acetic acid solution.
16. After 5-6 h of destaining, wash twice with distilled water and allow the gel (in water) to cool for 20 min in the freezer. This treatment improves the band quality as the gel shrinks a little at cold temperatures.

Photography:

17. Transfer the gel onto a glass plate, apply light from underneath and photograph.
18. Store the gel in Ziplock freezer bags or dry using commercial gel dryer.

Reference:

Hussain, A., H. Ramirez, W. Bushuk and W.M. Roca, 1986.

Field bean (*Phaseolus vulgaris* L.) cultivar identification by electrophoregrams of cotyledon storage proteins.
Euphytica, 35:729-732.

* Coomassie Blue Solution:

Dissolve 1g Coomassie Blue R-250 in 100 ml ethanol. Remove 4 ml of the above solution and add to 100 ml of Trichloro acetic acid. Mix and filter.

Cassava (Manihot esculenta Crantz)

Procedure

Preparation of material:

1. Cut stem portions from mature (11 months) plants and plant in sand : soil (1:1) mixture in greenhouse.
2. After 1 month, remove plants and wash roots with tap water.
3. Use root tip as source material for extraction of marker proteins.

Preparation of sample:

4. Weigh 0.5g of root tips.
5. Grind in porcelain mortar and pestle with 1 ml of the extraction buffer.*
6. Centrifuge at 27,000 Xg for 15 min at 10°C.
7. Use 20 uL of the clear supernatant for electrophoresis.

Preparation of gel:

8. Prepare 10% separating gel with 4% stacking as described in CIAT working document No.19 (1986) or as in LKB catalogue.

Electrophoresis:

Apparatus:	Bio Rad Protein II
Acrylamide concentration:	10%
Running buffer:	Tris borate, pH 9.0
Running voltage:	250 volts
Running time:	5 h
Running temperature	10°C

Staining:

9. Place gel into a stainless steel tray containing 200 ml of staining buffer** and incubate at 37°C for 10 min.

* Extraction buffer: Tris HCl 0.05M, pH 8.3 containing 5% PVP-40, 20% Sucrose, 0.5% Triton X 100 and 14 mM 2-mercaptoethanol.

** Staining buffer: Sodium phosphate 0.15M, pH 7.0.

10. Dissolve 100 mg α -naphthyl acetate and 50 mg β -naphthyl acetate in 5 ml acetone. Add little water till it is a bit turbid and pour it over the gel.
11. Incubate for 10 min at 37 °C.
12. Dissolve 150 mg of Fast Blue RR salt in small volume of water and filter. Pour this filtrate over the gel already floating in buffer and substrate solution.
13. Incubate at 37 C for further 1h until the bands become visible.
14. Stop the reaction by a rinse with water and fix with destaining solution.*

References:

Hussain, A. W. Bushuk, H. Ramirez and W.M. Roca, 1987.

Identification of cassava (*Manihot esculenta* Crantz) cultivars by electrophoretic patterns of esterase isozymes.
Seed Sci. & Technol. 15: 19-22.

Ramirez, H., A. Hussain, W. Roca and W. Bushuk, 1987.

Isozyme electrophoregrams of sixteen enzymes in five tissues of cassava (*Manihot esculenta* Crantz) varieties.
Euphytica 36: 39-48.

-
- * Destaining solution:
Methanol: Acetic acid : Water (6:1:14, v/v/v).

Pasture Legumes
(Centrosema macrocarpum, C. pubescens and C. acutifolium)

Procedure

Preparation of sample:

1. Remove husks manually or mechanically.
2. Grind seeds to a fine meal in mortar and pestle.
3. Weigh 0.1g of seed meal and add 1 ml of extraction solution* and vortex for 5 min.
4. Incubate at 40 C for 2h.
5. Centrifuge at 27,000 Xg for 20 min at 10 °C.
6. Transfer supernatant carefully into a vial and use 6 uL for electrophoresis.

Preparation of gel:

Preparation of acid gel is similar as described for field beans. Only Sucrose contents were increased from 5 to 10 g per gel and Hydrogen peroxide was reduced from 175 uL to 150 uL.

Electrophoresis:

Apparatus:	Bushuk's vertical unit
Acrylamide concentration:	10%
Running buffer:	Aluminum lactate, pH 3.1
Running current:	20 mA/gel
Running temperature:	20 °C

Staining/destaining and photography:

Same as described for field beans (see above).

Reference:

Hussain, A., H. Ramirez, W.M. Roca and W. Bushuk

Identification of cultivars of pasture legumes (Centrosema macrocarpum, C. pubescens and C. sp.n.) by acid polyacrylamide gel electrophoresis of cotyledon storage proteins.
Euphytica:

-
- * Extraction solution:
- | | |
|--|--------|
| a. Sodium acetate 0.5M, pH 5.8: | 1 part |
| b. Methyl green dye 0.5% in 40%
Sucrose solution: | 1 part |

Pasture Legume (Desmodium ovalifolium)

Procedure

Preparation of sample:

1. Grind seeds alongwith dry ice or using precooled mortar and pestle.
2. Weigh 0.4g of meal into a centrifuge tube.
3. Add 2 ml of 1% Sodium dodecyl sulphate solution.
4. Vortex for 5 min at room temperature.
5. Centrifuge at 27,000 Xg for 20 min at 5 °C.
6. Discard residue and heat supernatant in boiling water for 3 min.
7. Cool and centrifuge at 27,000Xg for 20 min at 5 °C.
8. Discard residue and add Aluminum lactate buffer, pH 3.1 to the supernatant in equal volume (1:1 v/v).
9. Mix and centrifuge at 27,000 Xg, 20 min at 5 °C.
10. Discard supernatant and wash residue with distilled water.
11. Redissolve residue in 0.4 ml of 1% SDS.
12. Centrifuge at 27,000 Xg for 20 min at 5 °C.
13. Decant supernatant, discard residue and add Sucrose to 10% level.
14. Add a drop of bromophenol blue and use 6 uL for electrophoresis.

Preparation of gel:

Same as described in LKB-manual or in CIAT Working Document No. 19. Gel does not contain SDS.

Electrophoresis:

Apparatus:	LKB 2001
Acrylamide concentration:	15%
Running buffer	Tris borate, pH 8.9
Running voltage	250 volts
Running time	7h
Running temperature	10 °C

Staining:

After electrophoresis stain gel with Coomassie Blue R-250 and destain with Methanol: Acetic acid: Water (6:1:14, v/v/v).

Reference:

Hussain, A., H. Ramirez, W. Bushuk and W.M. Roca.

Identification of Desmodium cultivars by electrophoresis.
Can. J. Plant Sci.

Pasture Legume (Stylosanthes capitata)

Procedure

Preparation of sample:

1. Remove seed coat by rubbing seeds against sand paper.
2. Crush 0.1g of seeds in 3 ml of the extraction buffer* using mortar and pestle.
3. Centrifuge at 12,000 Xg for 15 min at 10° C.
4. Use 15 uL of supernatant for electrophoresis.

Electrophoresis

Apparatus	LKB 2001
Acylamide concentration:	12%
Running buffer:	Tris borate, pH 9.0 (0.175M Tris + 0.0175M Boric acid)
Running voltage	250 volts
Running time	6h
Running temperature	8° C

Staining and photography:

Stain gel overnight with Coomassie Blue R-250 and destain with Methanol:Acetic acid: Water (6:1:14, v/v/v). Using diffused light from underneath , take photographs. Gel can be stored in ziplock freezer bags if necessary.

Reference:

Hussain, A., H. Ramirez, W.M. Roca and W. Bushuk,

Identification of cultivars of *Stylosanthes capitata* by polyacrylamide gel electrophoresis of seed proteins. *Euphytica*:

* Extraction buffer:

Tris HCl, pH 8.3 with 20% Sucrose and 0.002% Dithiothreitol.

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A. HUSSAIN, W. BUSHUK

Food Science Department,
University of Manitoba,
Winnipeg, Canada, R3T 2N2
(Publication No. 115)

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Referencias

-Hussain, A., W. Bushuk, H. Ramirez and W.M. Roca, 1987

Identification of cassava (Manihot esculenta Crantz) cultivars by electrophoretic patterns of esterase isozymes.
Seed Sci. & Technol. 15: 19-22

-Ramirez, H., A. Hussain, W. Roca and W. Bushuk, 1987

Isozyme electrophoregrams of sixteen enzymes in five tissues of cassava (Manihot esculenta Crantz) varieties.
Euphytica 36: 39-48

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Identification of cultivars of pasture legumes
(Centrosema macrocarpum, C. pubescens and C. sp.n)
by acid polyacrylamide gel electrophoresis of
cotyledon storage proteins.
Euphytica:

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Identification of Desmodium cultivars by
electrophoresis.
Can. J. Plant Sci.

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Referencia

Hussain, A., H. Ramirez, W.M. Roca and W. Bushuk.

Identification of cultivars of Stylosanthes capitata by polyacrylamide gel electrophoresis of seed proteins.
Euphytica:

- * Extraction buffer:
Tris HCl, pH 8.3, con 20% de sucrosa y 0.002% de Ditiotreitól.
- ** Buffer de corrida:
Trisma Base 0.175 M + Acido Borico 0.0175 M.